



SMZ161 Series

# SMZ-161-TL





SMZ161 Series

## SMZ-161-TL

<b>Optical system</b>	Greenough
<b>Observation tube</b>	Trinocular head
<b>Inclination</b>	45° inclined
<b>Trinocular light split</b>	Fixed 50:50
<b>Interpupillary distance</b>	50-75mm
<b>Diopter adjustment</b>	On both tubes, +/- 5 diopter
<b>Eyepieces</b>	Widefield WF10X/20mm
<b>Objectives system</b>	Zoom, ratio 6:1
<b>Objectives</b>	0.75X~4.5X
<b>Working distance</b>	110mm
<b>Stand type</b>	Pole type stand with transmitted illumination
<b>Head holder</b>	For Ø25mm column and Ø76mm head with incident illumination
<b>Focus mechanism</b>	Coarse focusing system with tension adjustment
<b>Focusing stroke</b>	43mm
<b>Incident illumination</b>	Quartz halogen 12V/10W with intensity control
<b>Transmitted illumination</b>	Quartz halogen 12V/20W with intensity control
<b>Transformer</b>	Internal
<b>Power supply</b>	100-240V (CE)
<b>Accessories included</b>	Dust cover, black/white and frosted stage plates and blue filter
<b>Dimensions LxWxH</b>	237x170x398mm
<b>Base</b>	237x170mm
<b>Column/Arm</b>	Height 205mm, Ø25mm
<b>Net weight</b>	3.7kg