

# US- $\alpha$ (optional accessory)



Prevent build-up on the prism surface with the ultrasonic cleaner.

\* The unit prevents sample build-up, it does not remove sample build-up after it has accumulated.


The US- $\alpha$  prevents build-up on the prism surface of the PRM-100 $\alpha$  and PRM-2000 $\alpha$ .  
As with other fittings, it can be affixed in the piping system.

## Specifications

Ultrasonic Cleaner US- $\alpha$

Cat.No.	For 1S fitting : 9112-1 For 2S fitting : 9112-2 For 3S fitting : 9112-3
Power supply	AC100 to 240V, 50/60Hz
Power consumption	5VA
Cable connector	Between US- $\alpha$ 's ultrasonic device part and power supply unit: (5m)
Dimensions and weight	Ultrasonic device part: 6.6×7.9×6.6cm, 1.4kg or more Power supply unit: 19.2×10×24cm, 3.5kg

## Connection example: US- $\alpha$ and PRM-100 $\alpha$



## Dimensions (unit of length: mm)

