

Thermo Scientific Thermo-Flask Benchtop Liquid Nitrogen Containers

Our Thermo-Flask liquid nitrogen containers feature inner vessels of borosilicate glass evacuated to 1×10^{-5} mm Hg. Available in stainless steel, and featuring vented clamped lids, carrying handles and rubber-cushioned bases.

- Durable, corrosion-resistant stainless steel exterior
- Shallow wide-mouth ideal for applications requiring flash freezing*
- Convenient dish-shaped design perfect for benchtop use
- “Push-fit” insulated lid, vented to prevent pressure build-up

> Thermo Scientific Benchtop Liquid Nitrogen Containers

Model No.	LN ₂ Capacity (liters)	Inside Neck Diameter in. (cm)	Exterior Dimensions Diameter x Height in. (cm)	Shipping Weight
2129	1	5.1 (13)	6.3 x 4.5 (16 x 11.4)	8 lbs. / 3.6 kg
2130	1.9	6.1 (15.5)	7.3 x 5.4 (18.5 x 13.7)	10 lbs. / 4.5 kg
2122	1.01	3.4 (8.6)	4.6 x 9 (11.7 x 22.9)	6 lbs. / 2.7 kg
2123	2.01	4.2 (10.7)	5.6 x 10.6 (14.2 x 26.9)	8 lbs. / 3.6 kg
2124	4.51	5.9 (15)	7.2 x 13.8 (18.3 x 35.1)	10 lbs. / 4.5 kg

> Accessory

Model No.	Description
AY509X6	12 mL Dipper



*The shallow wide mouth flasks do not have lids. Only the flasks with handles have lids.