



K Series

DSK700





K Series

DSK700

Optical system	Common Main Objective (CMO) Infinity optical system
Observation tube	Dual view Binocular head
Inclination	45° inclined
Interpupillary distance	54-76mm
Diopter adjustment	On both tubes, +/- 5 diopter
Eyepieces	Widefield WF10X/23mm
Objectives system	Zoom, ratio 5.2:1
Objectives	0.6X~3.1X
Working distance	89mm
Stand type	Pole type stand
Head holder	Optical bridge with light pointer
Focus mechanism	Coarse focusing system with tension adjustment
Transformer	External
Power supply	220-240V (CE)
Accessories included	Dust cover